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<b>GENERAL</b>		TEST RANGE LOW	5.8
MEASUREMENT MODE	ABSORB	HIGH	120
REACTION MODE	R-S <1>	NORM.RANGE LOW	13.3
CALIBRATION MODE	SLOPE AVG <2>	HIGH	63.9
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	END POINT<1>
WAVELENGTH	340nm<2>	READINGS FIRST	7 LAST : 20
DECIMAL POSITION	0	REACTION LIMIT	NO
UNIT	U / L <21>	<b>CALIBRATION</b>	
<b>ANALYSIS</b>		CALIB . INTERVAL	ON REQUEST <3>
POST DIL . FACTOR	NO <SPACE>	BLANK	
CONC . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
SAMPLE CYCLE	1	HIGH	NO <SPACE>
VOLUME	30	BLANK RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	30	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	300	STD - 2	NO
<b>CALCULATION</b>		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	DECREASE <2>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>ACE</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		CONVERS . FACTOR	1.00000
MEASUREMENT MODE	ABSORB	OFFSET	0.00000
REACTION MODE	R-S-SR1 <3>	TEST RANGE LOW	0
CALIBRATION MODE	SLOPE AVG <2>	HIGH	200
REAGENT BLANK	REAG / DIL <2>	NORM.RANGE LOW	0
CLEANER	NO <1>	HIGH	15
WAVELENGTH	550nm<4>	NUMBER OF STEPS	1
DECIMAL POSITION	1	CALC . STEP A	KINSEARCH
UNIT	U / L <21>	READINGS FIRST	12 LAST : 20
<b>ANALYSIS</b>		REACTION LIMIT	NO
POST DIL . FACTOR	NO <SPACE>	<b>CALIBRATION</b>	
CONC . FACTOR	NO <SPACE>	CALIB . INTERVAL	ON REQUEST <3>
SAMPLE CYCLE	1	BLANK	
VOLUME	6	REAG . RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	5.0	BLANK RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	200	STANDARD POS	#
START R1 CYCLE	8	STD - 1	*
VOLUME	100	STD - 2	NO
DILUTION NAME	H2O <00>	STD - 3	NO
VOLUME	5.0		
<b>CALCULATION</b>		REPLICATE	DUP <2>
SAMPLE LIMIT	NO <SPACE>	DEVIATION	10.0%
REAC . DIRECTION	INCREASE <1>	<b>ADA</b>	
CHECK	OFF		

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GENERAL		TEST RANGE LOW	1
MEASUREMENT MODE	ABSORB	HIGH	6
REACTION MODE	R-S <1>	NORM.RANGE LOW	3.5
CALIBRATION MODE	SLOPE AVG <2>	HIGH	5
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	END POINT <1>
WAVELENGTH	600nm <5>	READINGS FIRST	CB LAST : 10
DECIMAL POSITION	2	CALIBRATION	
UNIT	G / DL <11>	CALIB . INTERVAL	ON REQUEST <3>
ANALYSIS		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	3.0	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>	STANDARD POS	#
VOLUME	25.0	STD - 1	*
REAGENT CYCLE	1	STD - 2	NO
VOLUME	300	STD - 3	NO
CALCULATION			
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>ALBUMIN</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		CONVERS . FACTOR	1.00000
MEASUREMENT MODE	ABSORB	OFFSET	0.00000
REACTION MODE	R-S-SR1 <3>	TEST RANGE LOW	5
CALIBRATION MODE	SLOPE AVG <2>	HIGH	1000
REAGENT BLANK	REAG / DIL <2>	NORM.RANGE LOW	64
CLEANER	NO <1>	HIGH	306
WAVELENGTH	405nm <2>	NUMBER OF STEPS	1
DECIMAL POSITION	0	CALC . STEP A	KINSEARCH<3>
UNIT	U / L <21>	READINGS FIRST	8    LAST : 15
<b>ANALYSIS</b>		REACTION LIMIT	NO
POST DIL . FACTOR	NO <SPACE>	<b>CALIBRATION</b>	
CONC . FACTOR	NO <SPACE>	CALIB . INTERVAL	ON REQUEST <3>
SAMPLE CYCLE	1	BLANK	
VOLUME	5	REAG . RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	10.0	BLANK RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	240	STANDARD POS	#
START R1 CYCLE	5	STD - 1	*
VOLUME	60	STD - 2	NO
DILUTION NAME	H2O <00>	STD - 3	NO
VOLUME	5.0		
<b>CALCULATION</b>		REPLICATE	DUP <2>
SAMPLE LIMIT	NO <SPACE>	DEVIATION	10.0%
REAC . DIRECTION	INCREASE <1>	<b>ALP</b>	
CHECK	ON <1>		

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<b>GENERAL</b>		TEST RANGE LOW	5
MEASUREMENT MODE	ABSORB	HIGH	1000
REACTION MODE	R-S <1>	NORM.RANGE LOW	64
CALIBRATION MODE	SLOPE AVG <2>	HIGH	306
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINESEARCH <3>
WAVELENGTH	405nm <2>	READINGS FIRST	4 LAST : 11
DECIMAL POSITION	0	<b>CALIBRATION</b>	
UNIT	U / L <21>	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	5	HIGH	NO <SPACE>
DILUTION NAME	H2O <NO>	STANDARD POS	#
VOLUME	10	STD - 1	*
REAGENT CYCLE	1	STD - 2	NO
VOLUME	300	STD - 3	NO
<b>CALCULATION</b>			
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	ON <1>	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>ALP (WR) 4+1</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		CONVERS . FACTOR	1.00000
MEASUREMENT MODE	ABSORB	OFFSET	0.00000
REACTION MODE	R-S-SR1 <3>	TEST RANGE LOW	2
CALIBRATION MODE	SLOPE AVG <2>	HIGH	300
REAGENT BLANK	REAG / DIL <2>	NORM.RANGE LOW	0
CLEANER	NO <1>	HIGH	40
WAVELENGTH	340nm <1>	NUMBER OF STEPS	1
DECIMAL POSITION	1	CALC . STEP A	KINSEARCH
UNIT	U / L <21>	READINGS FIRST	8 LAST : 15
<b>ANALYSIS</b>		REACTION LIMIT	NO
POST DIL . FACTOR	NO <SPACE>	<b>CALIBRATION</b>	
CONC . FACTOR	NO <SPACE>	CALIB . INTERVAL	ON REQUEST <3>
SAMPLE CYCLE	1	BLANK	
VOLUME	30	REAG . RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	5.0	BLANK RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	240	STANDARD POS	#
START R1 CYCLE	5	STD - 1	*
VOLUME	60	STD - 2	NO
DILUTION NAME	H2O <00>	STD - 3	NO
VOLUME	5.0		
<b>CALCULATION</b>		REPLICATE	DUP <2>
SAMPLE LIMIT	NO <SPACE>	DEVIATION	10.0%
REAC . DIRECTION	DECREASE <2>	<b>ALT</b>	
CHECK	OFF		

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GENERAL		TEST RANGE LOW	2
MEASUREMENT MODE	ABSORB	HIGH	300
REACTION MODE	R-S <1>	NORM.RANGE LOW	0
CALIBRATION MODE	SLOPE AVG <2>	HIGH	40
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINSEARCH
WAVELENGTH	340nm <1>	READINGS FIRST	4 LAST : 11
DECIMAL POSITION	1	CALIBRATION	
UNIT	U / L <21>	CALIB . INTERVAL	ON REQUEST <3>
ANALYSIS		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	30	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>	STANDARD POS	#
VOLUME	30	STD - 1	*
REAGENT CYCLE	1	STD - 2	NO
VOLUME	300	STD - 3	NO
CALCULATION			
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	DECREASE <2>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>ALT (WR) 4+1</b>	
OFFSET	0.00000		

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GENERAL		TEST RANGE LOW	5
MEASUREMENT MODE	ABSORB	HIGH	2000
REACTION MODE	R-S <1>	NORM.RANGE LOW	0
CALIBRATION MODE	SLOPE AVG <2>	HIGH	100
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINESEARCH <3>
WAVELENGTH	405nm <2>	READINGS FIRST	3 LAST : 9
DECIMAL POSITION	0	CALIBRATION	
UNIT	U / L <21>	CALIB . INTERVAL	ON REQUEST <3>
ANALYSIS		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	5.0	HIGH	NO <SPACE>
DILUTION NAME	H2O <NO>	STANDARD POS	#
VOLUME	5.0	STD - 1	*
REAGENT CYCLE	1	STD - 2	NO
VOLUME	250	STD - 3	NO
CALCULATION			
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>a-AMYLASE</b>	
OFFSET	0.00000		

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GENERAL		REAC . DIRECTION	INCREASE <1>
MEASUREMENT MODE	ABSORB	CHECK	OFF
REACTION MODE	R-S-SR1 <3>	CONVERS . FACTOR	1.00000
CALIBRATION MODE	SLOPE AVG <2>	OFFSET	0.00000
REAGENT BLANK	REAG / DIL <2>	TEST RANGE LOW	1.9
CLEANER	NO <1>	HIGH	318
WAVELENGTH	340nm <1>	NORM.RANGE LOW	0
DECIMAL POSITION	0	HIGH	24
UNIT	U / L <21>	NUMBER OF STEPS	1
<b>ANALYSIS</b>		CALC . STEP A	KINSEARCH
POST DIL . FACTOR	NO <SPACE>	READINGS FIRST	9 LAST : 18
CONC . FACTOR	NO <SPACE>	REACTION LIMIT	NO
SAMPLE CYCLE	1	<b>CALIBRATION</b>	
VOLUME	24	CALIB . INTERVAL	ON REQUEST <3>
DILUTION NAME	H2O	BLANK	
VOLUME	5.0	REAG . RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	240	BLANK RANGE LOW	NO <SPACE>
START R1 CYCLE	4	HIGH	NO <SPACE>
VOLUME	60		
DILUTION NAME	H2O <00>	REPLICATE	DUP <2>
VOLUME	5.0	DEVIATION	10.0%
<b>CALCULATION</b>		<b>CK MB</b>	
SAMPLE LIMIT	NO <SPACE>		

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<b>GENERAL</b>		TEST RANGE LOW	1.9
MEASUREMENT MODE	ABSORB	HIGH	318
REACTION MODE	R-S <1>	NORM.RANGE LOW	0
CALIBRATION MODE	SLOPE AVG <2>	HIGH	24
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINSEARCH
WAVELENGTH	340nm <1>	READINGS FIRST	6 LAST : 15
DECIMAL POSITION	0	<b>CALIBRATION</b>	
UNIT	U / L <21>	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	24	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>	STANDARD POS	#
VOLUME	5.0	STD - 1	*
REAGENT CYCLE	1	STD - 2	NO
VOLUME	300	STD - 3	NO
<b>CALCULATION</b>			
SAMPLE LIMIT	NO <SPACE>	REPLICATE	DUP <2>
REAC . DIRECTION	INCREASE <1>	DEVIATION	10.0%
CHECK	OFF		
CONVERS . FACTOR	1.00000	<b>CK MB (WR) 4+1</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		CONVERS . FACTOR	1.00000
MEASUREMENT MODE	ABSORB	OFFSET	0.00000
REACTION MODE	R-S-SR1 <3>	TEST RANGE LOW	0.2
CALIBRATION MODE	SLOPE AVG <2>	HIGH	15
REAGENT BLANK	REAG / DIL <2>	NORM.RANGE LOW	0.6
CLEANER	NO <1>	HIGH	1.4
WAVELENGTH	505nm <5>	NUMBER OF STEPS	1
DECIMAL POSITION	2	CALC . STEP A	KINSEARCH
UNIT	MG / DL	READINGS FIRST	4 LAST : 7
<b>ANALYSIS</b>		REACTION LIMIT	NO
POST DIL . FACTOR	NO <SPACE>	<b>CALIBRATION</b>	
CONC . FACTOR	NO <SPACE>	CALIB . INTERVAL	ON REQUEST <3>
SAMPLE CYCLE	1	BLANK	
VOLUME	20	REAG . RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	10	BLANK RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	100	STANDARD POS	#
START R1 CYCLE	2	STD - 1	*
VOLUME	90	STD - 2	NO
DILUTION NAME	H2O <00>	STD - 3	NO
VOLUME	5.0		
<b>CALCULATION</b>		REPLICATE	DUP <2>
SAMPLE LIMIT	NO <SPACE>	DEVIATION	10.0%
REAC . DIRECTION	INCREASE <1>	<b>CREATININE JAFFE</b>	
CHECK	OFF		

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<b>GENERAL</b>		TEST RANGE LOW	0.2
MEASUREMENT MODE	ABSORB	HIGH	15
REACTION MODE	R-S <1>	NORM.RANGE LOW	0.6
CALIBRATION MODE	SLOPE AVG <2>	HIGH	1.4
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINSEARCH
WAVELENGTH	505nm <3>	READINGS FIRST	2 LAST : 5
DECIMAL POSITION	2	<b>CALIBRATION</b>	
UNIT	MG / DL	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	25	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	10	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	250	STD - 2	NO
<b>CALCULATION</b>		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>CREATININE JAFFE (WR) 1+1</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		CONVERS . FACTOR	1.00000
MEASUREMENT MODE	ABSORB	OFFSET	0.00000
REACTION MODE	R-S-SR1 <3>	TEST RANGE LOW	5
CALIBRATION MODE	SLOPE AVG <2>	HIGH	300
REAGENT BLANK	REAG / DIL <2>	NORM.RANGE LOW	7
CLEANER	NO <1>	HIGH	50
WAVELENGTH	405nm <2>	NUMBER OF STEPS	1
DECIMAL POSITION	1	CALC . STEP A	KINSEARCH
UNIT	U / L <21>	READINGS FIRST	9 LAST : 18
<b>ANALYSIS</b>		REACTION LIMIT	NO
POST DIL . FACTOR	NO <SPACE>	<b>CALIBRATION</b>	
CONC . FACTOR	NO <SPACE>	CALIB . INTERVAL	ON REQUEST <3>
SAMPLE CYCLE	2	BLANK	
VOLUME	30	REAG . RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	10	BLANK RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	240	STANDARD POS	#
START R1 CYCLE	4	STD - 1	*
VOLUME	60	STD - 2	NO
DILUTION NAME	H2O <00>	STD - 3	NO
VOLUME	5.0		
<b>CALCULATION</b>		REPLICATE	DUP <2>
SAMPLE LIMIT	NO <SPACE>	DEVIATION	10.0%
REAC . DIRECTION	INCREASE <1>	<b>g GT</b>	
CHECK	OFF		

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GENERAL		TEST RANGE LOW	5
MEASUREMENT MODE	ABSORB	HIGH	300
REACTION MODE	R-S <1>	NORM.RANGE LOW	7
CALIBRATION MODE	SLOPE AVG <2>	HIGH	50
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINSEARCH
WAVELENGTH	405nm <2>	READINGS FIRST	6 LAST : 15
DECIMAL POSITION	1	CALIBRATION	
UNIT	U / L <21>	CALIB . INTERVAL	ON REQUEST <3>
ANALYSIS		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	25	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	10.0	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	250	STD - 2	NO
CALCULATION		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>g GT (WR) 4+1</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		TEST RANGE LOW	0
MEASUREMENT MODE	ABSORB	HIGH	140
REACTION MODE	R-S <1>	NORM.RANGE LOW	4.5
CALIBRATION MODE	SLOPE AVG <2>	HIGH	19.8
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	END POINT <1>
WAVELENGTH	505nm <3>	READINGS FIRST	CB LAST : 13
DECIMAL POSITION	0	<b>CALIBRATION</b>	
UNIT	MG / DL <12>	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	3.0	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	10	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	300	STD - 2	NO
<b>CALCULATION</b>		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>LACTATE</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		CONVERS . FACTOR	1.00000
MEASUREMENT MODE	ABSORB	OFFSET	0.00000
REACTION MODE	R-S-SR1 <3>	TEST RANGE LOW	10
CALIBRATION MODE	SLOPE AVG <2>	HIGH	1600
REAGENT BLANK	REAG / DIL <2>	NORM.RANGE LOW	0
CLEANER	NO <1>	HIGH	530
WAVELENGTH	340nm <1>	NUMBER OF STEPS	1
DECIMAL POSITION	0	CALC . STEP A	KINSEARCH
UNIT	U / L <21>	READINGS FIRST	8 LAST : 13
<b>ANALYSIS</b>		REACTION LIMIT	NO
POST DIL . FACTOR	NO <SPACE>	<b>CALIBRATION</b>	
CONC . FACTOR	NO <SPACE>	CALIB . INTERVAL	ON REQUEST <3>
SAMPLE CYCLE	1	BLANK	
VOLUME	6	REAG . RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	10	BLANK RANGE LOW	NO <SPACE>
REAGENT CYCLE	1	HIGH	NO <SPACE>
VOLUME	240		
START R1 CYCLE	5	STANDARD POS	#
VOLUME	60	STD - 1	*
DILUTION NAME	H2O <00>	STD - 2	NO
VOLUME	5.0	STD - 3	NO
<b>CALCULATION</b>		REPLICATE	DUP <2>
SAMPLE LIMIT	NO <SPACE>	DEVIATION	10.0%
REAC . DIRECTION	DECREASE <2>	<b>LDH</b>	
CHECK	OFF		

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<b>GENERAL</b>		TEST RANGE LOW	10
MEASUREMENT MODE	ABSORB	HIGH	1600
REACTION MODE	R-S <1>	NORM.RANGE LOW	0
CALIBRATION MODE	SLOPE AVG <2>	HIGH	530
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINSEARCH
WAVELENGTH	340nm <1>	READINGS FIRST	4 LAST : 10
DECIMAL POSITION	0	REACTION LIMIT	NO
UNIT	U / L <21>	<b>CALIBRATION</b>	
<b>ANALYSIS</b>		CALIB . INTERVAL	ON REQUEST <3>
POST DIL . FACTOR	NO <SPACE>	BLANK	
CONC . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
SAMPLE CYCLE	1	HIGH	NO <SPACE>
VOLUME	6	BLANK RANGE LOW	NO <SPACE>
DILUTION NAME	H2O <00>	HIGH	NO <SPACE>
VOLUME	10		
REAGENT CYCLE	1	STANDARD POS	#
VOLUME	300	STD - 1	*
<b>CALCULATION</b>		STD - 2	NO
SAMPLE LIMIT	NO <SPACE>	STD - 3	NO
REAC . DIRECTION	DECREASE <2>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>LDH 4+1</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		TEST RANGE LOW	0.5
MEASUREMENT MODE	ABSORB	HIGH	5
REACTION MODE	R-S <1>	NORM.RANGE LOW	1.8
CALIBRATION MODE	SLOPE AVG <2>	HIGH	2.6
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	END POINT <1>
WAVELENGTH	550nm<4>	READINGS FIRST	CB LAST :10
DECIMAL POSITION	2	<b>CALIBRATION</b>	
UNIT	MG / DL <11>	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	3	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	25.0	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	300	STD - 2	NO
<b>CALCULATION</b>		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	ON <1>	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>MAGNESIUM</b>	
OFFSET	0.00000		

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GENERAL		TEST RANGE LOW	1
MEASUREMENT MODE	ABSORB	HIGH	14
REACTION MODE	R-S <1>	NORM.RANGE LOW	6.6
CALIBRATION MODE	SLOPE AVG <2>	HIGH	8.8
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	END POINT <1>
WAVELENGTH	550nm <3>	READINGS FIRST	CB LAST 10
DECIMAL POSITION	0	<b>CALIBRATION</b>	
UNIT	G / DL <12>	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	7	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	10	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	280	STD - 2	NO
<b>CALCULATION</b>		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>TOTAL PROTEIN</b>	
OFFSET	0.00000		

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<b>GENERAL</b>		TEST RANGE LOW	5
MEASUREMENT MODE	ABSORB	HIGH	700
REACTION MODE	R-S <1>	NORM.RANGE LOW	0
CALIBRATION MODE	SLOPE AVG <2>	HIGH	200
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	END POINT <1>
WAVELENGTH	505nm <3>	READINGS FIRST	CB LAST 20
DECIMAL POSITION	0	<b>CALIBRATION</b>	
UNIT	MG / DL <12>	CALIB . INTERVAL	ON REQUEST <3>
<b>ANALYSIS</b>		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	3	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	10	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	300	STD - 2	NO
<b>CALCULATION</b>		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	INCREASE <1>	REPLICATE	DUP <2>
CHECK	OFF	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>TRIGLYCERIDES</b>	
OFFSET	0.00000		

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GENERAL		TEST RANGE LOW	2
MEASUREMENT MODE	ABSORB	HIGH	200
REACTION MODE	R-S <1>	NORM.RANGE LOW	15
CALIBRATION MODE	SLOPE AVG <2>	HIGH	44
REAGENT BLANK	REAG / DIL <2>	NUMBER OF STEPS	1
CLEANER	NO <1>	CALC . STEP A	KINSEARCH
WAVELENGTH	340nm <1>	READINGS FIRST	2 LAST : 5
DECIMAL POSITION	1	CALIBRATION	
UNIT	MG / DL <12>	CALIB . INTERVAL	ON REQUEST <3>
ANALYSIS		BLANK	
POST DIL . FACTOR	NO <SPACE>	REAG . RANGE LOW	NO <SPACE>
CONC . FACTOR	NO <SPACE>	HIGH	NO <SPACE>
SAMPLE CYCLE	1	BLANK RANGE LOW	NO <SPACE>
VOLUME	3	HIGH	NO <SPACE>
DILUTION NAME	H2O <00>		
VOLUME	10	STANDARD POS	#
REAGENT CYCLE	1	STD - 1	*
VOLUME	300	STD - 2	NO
CALCULATION		STD - 3	NO
SAMPLE LIMIT	NO <SPACE>		
REAC . DIRECTION	DECREASE <1>	REPLICATE	DUP <2>
CHECK	ON <1>	DEVIATION	10.0%
CONVERS . FACTOR	1.00000	<b>UREA (WR) 4+1</b>	
OFFSET	0.00000		

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